



Winpact Solid State Fermentation System, FS-V-SA05P

Solid state fermentation (SSF) can be used for enzyme, antibiotics, biofuel, and organic acid production in the food, pharmaceutical, cosmetic, industries, etc. One of the features for Solid state fermentation is to create low water level cultivating conditions for fungus, mold, filamentous fungi, and some bacteria growth.

Winpact Solid State Fermentation system is designed for the laboratory scale research to get excellent results. It is featured with a 10.4" color touch screen, user-friendly interface and 4 built-in peristaltic pumps on the Linux based operation system. An automatic vessel angle control mechanism provides an outstanding mixing efficiency for solid state material research.

This system is suitable for both aerobic and anaerobic fermentation with three kinds of impellers available (Broken, Anchor and Spiral type).

(FS-V-SA05P)



*0°-90° rotation,120° for harvest

Features

- Fully integrated system specifically designed for solid-state fermentation research involving saccharification, hydrolysis and more.
- Programmable angle-adjustable (0-90° for culture control, 120° for harvest control) vessel tiling and stirring mechanism permits superior sample homogeneity
- Impellers are designed to reduce stickiness and it ensures sample integrity during the fermentation process.
- Integrated motor shaft & air sparger assembly creates precise, disturbance-free controls of aeration and mixing
- Chemically resistant double jacketed borosilicate glass vessel design
- Can be used with pH and DO probes to control culture conditions (anchor type impellers only)
- Customizable impellers and aeration controller available
- · Winpact Humidifier is available for real-time monitoring and humidity adjustment.

^{**}The minimum speed varies from 1-5 rpm depending on the medium viscosity.



www.majorsci.com info@majorsci.com





Impeller Type:









Broken

Anchor

Spiral

(FS-O-HMD)

Control Unit	Control Panel	10.4" color touch-screen Interface, (Resolution: 800 x 600 pixels)
	Communication Port	Remote control through Ethernet, Analog AUX port for system extension
	Storage Program	Up to 59,994 programs for different kinds of condition.
	Data Internal Storage	Up to 100 data files.
	Data External Storage Interface	USB port
	Cabinet Material	Front panel: ABS / Housing: Painted iron
	Rated Voltage	110V~/ 220V~ ; 50/60 Hz
Aeration	Inlet Gas Flow-meter	0, 1 – 10 LPM
Dimension	Overall Diameter 35	Omm; Overall Height with Condenser 683 mm; Overall Height without Condenser 448 mm
	Dimension (with vessel holder) 430mm (L) x 730mm (W) x 780 mm (H)	
Temperature	Heating	Thermostat system: Built-in heat exchanger, 550W heater/water circulation pump
	Cooling	Automatic cooling water valve
	Range	5°C (41°F) above coolant up to 60°C (140°F)
	Resolution	0.1°C
	Control Mode	Manual or programmable 15-step PID control.
Agitation	Drive	Removable top brushless motor
	Speed Range	0, 1 – 60 rpm
	Resolution	1rpm
	Control Mode	Manual or programmable 15-step PID control.
	Impeller	1. Broken type (FS-A-IM305)
		2. Anchor type (FS-A-IM408)
		3. Spiral type (FS-A-IM507)
		(Select one from the above type, and only anchor impeller can be used with pH and DO probes)
		*Note: Customized impellers are available.
		**In pH and DO measurement condition, the minimum medium volume is 4L and water content is
		more than 50%, tilting angle not over 30 degree.
		***The measure value of pH and DO may not accurate when using in solid-state culture condition.
		****pH and DO probe is not within the scope of warranty when using in solid-state vessel.
Vessel Swing	Angle Range	Normal operation: 0°~90°, adjustable time interval
		Harvest mode: 0°/ 120°
	Control Mode	Programmable control

INNOVATIVE CULTIVATION SOLUTIONS



www.majorsci.com info@majorsci.com

- *For more information, please contact your local distributors.
- * All images are for reference only, actual products might differ from the pictures above.
- * Technical specifications subject to change without notice.

Taiwan Headquarters

Headquarters:

No. 156, Sec. 1, Guoji Rd., Taoyuan Dist., Taoyuan City 33061, Taiwan T/+886-3-3762878 F/+886-3-3761310

Factory: No.19, Ln. 207, Huakang St. Bade Dist., Taoyuan City 33464 Taiwan T/ +886-3-3623319 F/ +886-3-3623133

US Office

19959 Sea Gull Way Saratoga, CA 95070 U.S.A. T/+1-408-366-9866

F/ +1-408-446-1107

Shanghai Office

Room 612, International business exhibition center, 9300 Hunan Road, Pudong, Shanghai, China National toll-free No.:400-823-9177 T/ +86-21-50795277 F/ +86-21-50795277

India Office

D. No. 12-13-99, Satguru Apts. Extn. Street. No. 3, Lane No.1, Tarnaka Secunderabad – 500 017. India T/ +91-40-27001515 T/ +91-40-27001586